



AB GERMA®

VACUUM MATTRESS MODEL AS 192



PRODUCT SHEET

AS 192 Vacuum mattress is a Full Body Vacuum Mattress that is developed and improved with the original high quality standards from Germa materials and vacuum products.

When immobilisation is needed, the vacuum mattress is easy to prepare due to Germa unique inner chamber system to make it easy adjustable, ready for use immediately after short preparation for premium immobilisation.

The mattress is equipped with 8 carry handles that are designed for comfortable carrying of a patient. Four restraints across the mattress to secure patient on top of the mattress.

VACUUM MATTRESS MODEL AS 192



PRODUCT SPECIFICATIONS

Approved according to European Standard
EN 1865
MDR 2017/745 Class 1 Medical device.

PRODUCT INFO

- 4 handles on each side – carry capacity 200 Kg
- Improved chambers system
- Includes removable, contoured pillow

Warranty period: 2 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated,
- flame retardant
- Inner bag: Polyester
- Granules: Polystyrene granules (flame retardant)

SPECIFICATIONS

- Length: 210 cm
- Width: 80 cm
- Height: 5 cm
- Weight 9 kg
- Loading capacity: 200 kg

TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -30°C to +70°C

ORDER INFORMATION

Art no: 271404000

Vacuum Mattress Model AS 192

SPARE PARTS

Art no: 12633800000

Art: Germa vacuum valve Black

Art no: 15000000111

Art: Repair Kit Mattresses

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se