



**AB GERMA®**

## DENTACOMFORT VACUUM PILLOW



# PRODUCT SHEET

The DentaComfort pillow has been developed together with the input from practising dentist to get the best possible product when it comes to shape, firmness and function. The DentaComfort pillow is used by dentists for the positioning of patients in the dental chair. It helps patients relax and makes the dentists job easier as the patient can be secured in a specific position with great comfort for the neck and head. The DentaComfort pillow can also give support for the dentist arms during procedures.

## DENTACOMFORT VACUUM PILLOW



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- The Vacuum Pillows consist of an airtight shell containing granule filled inner bag.
- The DentaComfort is shaped as a triangle for best positioning and function.
- Comes in two colours, Blue and white

Warranty period: 2 years

### MATERIAL

- Outer material: PVC-foil, flame retardant
- Inner bag: Nylone
- Granules: Polystyrene

Width: 50 cm  
Depth: 48 cm  
Hight: 5 cm  
Weight: 0,8 kg

### ORDER INFORMATION

Art no: 22393200000  
Positioning Pillow DentaComfort White

Art no: 22393240000  
Positioning Pillow DentaComfort Blue

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)  
Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)