



**AB GERMA®**

EASY HANDLE SET



# PRODUCT SHEET

Specific designed tools to support ergonomics, safety and comfortable handling of the patient in need for transport. Studies and development of ergonomic lifting technic is developing and improving, important to take the opportunities to make use of proven recommendations and lifting technic offered today.

In partnership with ergonomic specialists focused on development of training and tools, for best practice for paramedics when handling patients, which today are heavier than before in average.

Most important for the paramedics is to make use of ergonomically optimal body positioning and procedures, which often require professional input and training. In combination with well-designed tools, you can reduce risk of back injuries significantly.

# EASY HANDLE SET



X-RAY MRI AND  
CT TRANSLUCENT

## PRODUCT SPECIFICATIONS

### PRODUCT INFO

- 4 handles in one set, packed in a bag

Warranty period: 1 years

### MATERIAL

- Straps: Polyester
- Handles: PVC
- Bag: Polyester

### SPECIFICATIONS

- Length: 42 cm
- Width: 2,5 cm
- Thickness Handles: 3,3 cm
- Weight: 0,45 kg

### ORDER INFORMATION

Art no: 24271000737

Easy Handle Set

Weight: 0,45 kg

## BENEFITS

Patient handling tool  
Improved ergonomic patient handling, lifting and/  
or sliding

Lifting from knee level—reducing risk for back injury  
for the rescuer by 60%

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

www.germa.se  
Phone: +46 (0)44 12 30 30  
info@germa.se