



AB GERMA®

## EASYFIX PLUS II MOD VACUUM MATTRESS



# PRODUCT SHEET

Germa innovative revolutionary designed line of Vacuum Mattresses for the Military: EasyFix PLUS II MOD is a development work based on earlier experiences from the original high quality materials traditionally used producing Germa world class vacuum mattresses.

EasyFix PLUS II MOD vacuum mattress using improved body shape design to optimize the vacuum mattress capabilities. Offers 20% saving on storage space of the mattress compared to our All In One Vacuum Mattresses. EasyFix PLUS II MOD unique construction which deliver improved fixation of upper back, neck when excluding cervical collar on the patient and reinforced base plate that cover all of cervical spine. This together with fixation restraint system using X-system to support and create the most effective immobilisation.

## EASYFIX PLUS II MOD VACUUM MATTRESS



REINFORCED WOOD  
INNER SYSTEM



FLAME RETARDANT



HF-VELDED SEAMS



X-RAY MRI AND  
CT TRANSLUCENT



MULTI FIT VALVE WITH LID

## PRODUCT SPECIFICATIONS

Approved according to European Standard  
EN 1865 and EN 1789  
MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- 5 handles on each side – carry capacity 250 Kg
- Non IR reflective PVC
- New improved chambers system
- Improved design of neck support to support patient without cervical collar
- Reinforced reverse of the mattress to reduce friction and improve possible handling
- Improve head strap system for better grip
- Light weighted and small storage

Warranty period: 3 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, Flame retardant and Non IR reflective
- Inner bag: Polyester
- Granules: Polystyrene granules (flame retardant)
- Wood board upper part: Birch plywood

### SPECIFICATIONS

- Loading capacity: 250 kg
- Length: 208 cm
- Width: 130/72 cm (Head part)
- Thickness: 4,5 cm
- Weight: 9,6 kg (excluding pump etc.)

### TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -30°C to +70°C

### ORDER INFORMATION

Art no: 271511181

Vacuum Mattress EasyFix PLUS II MOD

Ready2Go Kit

– including restraint system, footrest Vacuum Pump, repair kit, bag

Weight: 12,5 kg

### SPARE PARTS

Art no: 12633800000

Art: Germa vacuum valve Black

Art no: 15000000113

Art: EasyFix MOD Repair Kit

### ACCESSORIES

Art no: 15000000117

Art: Step by Step Transfer Black II

Art no: 15000000034

Art: Easy Handle Set

Visiting address

AB Germa

Industrig. 54-56

SE-291 36 Kristianstad, Sweden

Post address

AB Germa

Box 2127

SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)

Phone: +46 (0)44 12 30 30

[info@germa.se](mailto:info@germa.se)

## BENEFITS

Made to correspond to EN 1789 and 1865.  
MDR 2017/745 Class I Medical device

## EASYFIX PLUS II MOD VACUUM MATTRESS

**Germa head pillow design including rigid wood part solution.** Efficient, quick and improved neck and head immobilisation, with visible ears in one go.

Support upper back, neck and head for safe immobilization when not using cervical collar

One-piece back of the mattress, providing the patient with even spine support when lifted.

Design of the upper part allows the rescuer to place the arm(s) of the patient in two different positions. This enables the rescuer access to the patient's arm for treatment.

Reinforced reverse of the mattress to reduce friction.

Made from plastic makes it ideal for X-ray. Allows preliminary X-ray diagnostics, direct in to x-ray no extra handling of patient.

Germa improved chamber technology for quick and easy preparation.  
Give rescuers quick handling and excellent immobilisation for the patient ensuring the granulate is evenly dispersed to give the patient stable and optimal immobilisation.

Handles with ergonomic grip, making patient lifts as comfortable as possible for rescuers.

Handles fixed to the bottom of mattress on piece, maintain rigid immobilisation when lifting the patient.

Robust welding—fully welded without any open seams or textiles, making hygienic cleaning of the mattress easy.

Valve with special designed hole, ensuring most suction devices fit directly into the valve, no adapters are needed.

Reverse air lock system, no closing of valve at stressful site.

Improved V-shaped design to make sure upper wide part will make all sized patients the very best immobilisation of the back, neck and head.

X-designed restraint belt fixation close to the patient, support immobilisation process and comfortable for patient. Supplied with green colour straps with fix lock, helping the rescuer to place the straps correctly and easy to stretch and adjust.

Outer material is made of flame retardant PVC-coated fabric, non IR reflective.

PVC coated fabric makes it easy to clean and disinfect with most common disinfectants, easy to keep clean and ready for use.

Robust and efficient material to hold vacuum for more than 24 hours, safety for patient and make confiding immobilisation.

High load capacity 250 kg, safe to handle heavy duty patients.

Repair Kit easy to use for effective repair, quick to repair no return needed.

Folds easily, 20% less storage space required compared to most standard vacuum mattresses.

High quality footrest aluminium vacuum pump evacuate with more than -700 hpa, durable and quick and time saving.