



AB GERMA®

EVAC CHAIR BACK SUPPORT



PRODUCT SHEET

EVAC CHAIR BACK SUPPORT

The Evac Chair Back Support prevents the patients head to fall through in between the handle bar when extended. The Evac Chair Back Support gives the patient a safer ride, easier to relax and be comfortable.

The new Evac Chair Back Support is easily attached to the chair with hooks & loops straps and it flawlessly adjusts to the handles when extended or shortened. paramedics in their work.

PRODUCT SPECIFICATIONS

PRODUCT INFO

- Made from a strong and sturdy mesch material.
- Has 6 hooks & loops straps for secure fastening on the chair.
- Low weight
- Easy to clean
- Give a safer transport for the patient
- The Mesch should not be folded but should be rolled together.
- Clean with a mild detergent and whipe of with wet cloth

Warranty period: 1 years

MATERIAL

- Mesch: PVC coated Polyester
- Straps: Polyester

SPECIFICATIONS

- Length: 57 cm
- Width: 22/13 cm
- Weight: 0,1 kg

ORDER INFORMATION

Art no: 15000000126

Evac Chair Back Support

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se