



**AB GERMA®**

**STRETCHER MATTRESS  
F2 4P MD**



## PRODUCT SHEET

The F2 Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient on a Ferno F2 stretcher/cot.

This quality mattress is designed with help of trauma teams to secure efficient handling with safety and comfort for the patient.

## STRETCHER MATTRESS F2 4P



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

### PRODUCT SPECIFICATIONS

#### PRODUCT INFO

Designed to fit on top of Ferno F2 stretcher/cot, fixation with hooks and loops.

Warranty period: 2 years

#### MATERIAL

- Outer material: Padova PVC - Polyester fibre, PVC-coated, flame retardant.
- Inner material: 70 mm Polyurethane Foam.

#### SPECIFICATION FOR THE F2 MATTRESS

- Length: 191 cm
- Width: 49 cm
- Height: 8 cm
- Weight 3,3 kg

#### ORDER INFORMATION

Art no: 271374000

Stretcher Mattress F2 4p Black

### BENEFITS

MDR 2017/745 Class 1 Medical device

Easy to assembly and adjust stretcher restraint system according to standard regulation EN 1789.

Designed to fit on top of Ferno F2 stretcher/cot, fixation with hooks and loops.

The PVC fabric is made strong and durable for long and hard use.

Extra comfort by using a combination with Comfort foam which reads and support patients body and even out pressure to be able to deliver premium safe and comfort feeling when laying on the mattress.

Hygienic, all welded product, easy to clean with most common disinfectants.