



**AB GERMA®**

## MINI TRANSFER MATTRESS MOD



# PRODUCT SHEET

Germa have together with our customers developed a very handy, flexible and useful patient handling tool e.g. to move a patient from a stretcher to a bed or lift a patient from the ground or narrow area onto for example a basket stretcher.

The extra wide handle part makes it easy to lift the patient in or out from a basket stretcher without hurting any fingers.

## MINI TRANSFER MATTRESS MOD



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device.

### PRODUCT INFO

- 4 handles on each side, 2 at head end
- Glass fiber back side

Warranty period: 2 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant and Non IR reflective
- Inner material: 10 mm thick Polyurethane Foam

### SPECIFICATION FOR THE TRANSFER MATTRESS

- Length: 136 cm
- Width: 78 cm
- Height: 1,5 cm
- Weight 3,0 kg

### TEST SPECIFICATIONS

- Working temperature of the mattress is -30°C to +70°C

### ORDER INFORMATION

Art no: 271048000

MINI TRANSFER MATTRESS MOD

### SPARE PARTS

Art no: 15005002817

Art: Repair Kit

## BENEFITS

## MINI TRANSFER MATTRESS MOD

One big feature you will experience using Mini Transfer is that it fits great in a basket stretcher, the patient can then lay more comfortable and safer in the basket stretcher.

When patient has to be rescued from a narrow or difficult environment e.g. the forrest, the smaller size makes it easy to use and to manouver around.

The positions of the handles makes it possible to carry/move the patient with great flexibility.

To make usage as easy as possible, the back of the Mini Trans mattress is made with friction reducing material to slide easy underneath the patient.

Glass fibre material on the back is very durable and can be used sliding on floor or ground.