



PRODUCT SHEET

Germa Motion Restriction Pillow,
vacuum pillow for Stretchers.

MRP Vacuum Pillow is designed to support and keep the patient's head in a stable position during ambulance transport. Medical evaluation clearly show positive effects for the patient and the paramedic's work performed during ambulance transport. MRP Vacuum Pillow is very comfortable in a non evacuated state but can quick be evacuated to support the neck and head.

The Germa manufactured MRP Pillow is intended for use together with the Germa Strecher/Transfer Mattresses and it fits most common stretchers on the market.

The flexible pillow is filled with a mix of foam, pillow fibres and granulate chambers to allow the patient to feel very comfortable. The MRP Vacuum Pillow is attached to the strecher itself with two straps.

MOTION RESTRICTION PILLOW (MRP)



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- Light weighted and small storage
- Warranty period: 1 year

MATERIAL

- Outer material: Polyester fibre, PU coated (flame retardant)
- Granules: Styroporgranules (flame retardant)
- Foam and polyester pillow filling.

SPECIFICATIONS SIZE:

- High: 36 cm
- Width: 48 cm
- Thickness: 10 cm (not evacuated)
- WEIGHT: 0,6 kg

ORDER INFORMATION:

Art: Germa Motion Restriction Pillow, vacuum, Black
Art no: 15000000115

ACCESSORIES:

Art: One-hand-operated Vacuum Pump
Art no: 272000701
Weight: 0.6 kg

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

MOTION RESTRICTION PILLOW (MRP)

When there is a suspected spinal injury the MRP vacuum pillow with it's specific construction has the intended use to limit the motion of the patient upper spine - to support the patient head and neck during transport to the hospital.

The pillow can be used as a normal comfort pilllow but when neck/spine support is needed it easily transforms to a rigid support.

Straps fixed to the head of the strecher

Robust welding—fully welded without any open seams or textiles, making hygienic cleaning of the pillow easy.

Valve with special designed hole, ensuring almost all suction devices fit directly into the valve, no adapters are needed. Reverse air lock system, no closing of valve at stressful site. Made from plastic makes it ideal for X-ray.

PU coated fabric makes it easy to clean and disinfect with most common disinfectants, easy to keep clean and ready for use.

High quality onehand operated vacuum pump evacuate with more than -400 hpa, durable and quick and time saving. According to EN 1865.





AB GERMA®

MOTION RESTRICTION PILLOW (MRP)

