



AB GERMA®

SPINE BOARD & SCOOP STRETCHER PADDING SET



PRODUCT SHEET

The Padding Set contents of 5 separate parts, one to be placed between the patient legs and two on each long side of the patient. Padding Set is intended to be used for immobilisation of a patient on top of a Spine Board or Scoop stretcher.

The Padding Set is developed to improve actual immobilisation, support and comfort of the patient when placed on a Spine Board or Scoop Stretcher. By using the Padding Set on both side of the body and between the leg a secure positioning and stabilisation is achieved during transport. For superiour security the Xpider Strap is recommended to be used as well.

SPINE BOARD & SCOOP STRETCHER PADDING SET

PRODUCT SPECIFICATIONS

PRODUCT INFO

- Premium immobilisation and optimal positioning of the patient on top of a Spine Board/Scoop Board together with the Xpider Strap.
- Easy to adjust for optimal patient positioning, limit pressure points on the patient.
- Colour coded for side and leg support paddings
- Easy to use and designed for easy cleaning.

Warranty period: 2 years

MATERIAL

- Orange: PVC—polyester fabric material (flame retardant).
- Blue: PVC foil (flame retardant).
- Foam Polyurethane
- Granules: Polystyrene granules (flame retardant)
- Bag: Polyamid

SPECIFICATIONS

- Length: 70 cm
- Diameter: 9 cm
- Weight: 0,3 kg
- Length: 50 cm
- Hight: 8/8,5 cm
- Weight: 0,225 kg

Weight of padding set: 1,2 kg

ORDER INFORMATION

Art no: 22221800001
Padding 500x80x85mm

Art no: 22221800002
Padding roll 70cm diameter 9cm

Art no: 22221800010
Padding set (4+1)