



AB GERMA®

# RESCUE LITE FULL BODY SPLINT MOD



# PRODUCT SHEET

Intended use is to give good stabilization for the patient being rescued in combination with e.g. a basket stretcher. Light construction for rescuers to carry with them in difficult terrain.

Rescue Lite Full Body Splint MOD's design is based on feedback from medical and rescue experts with great experience from rescue operations in difficult terrain. Construction based on original high quality materials traditionally used producing Germa world class vacuum mattresses.

Rescue Lite Full Body Splint MOD is utilizing improved body shape design to optimize the vacuum full body splint capabilities. Offers 35% saving on storage space compared to Germa Easyfix Plus II MOD Vacuum Mattresse.

## RESCUE LITE FULL BODY SPLINT MOD

The Rescue Lite Full Body Splint features a unique design that provides effective stabilization of the entire body, specifically the upper back and neck, while excluding the cervical collar. It also includes a reinforced base plate that covers the entire cervical spine. This together with fixation restraint system using X-system to support and create the most effective immobilisation.



FLAME RETARDANT



X-RAY MRI AND  
CT TRANSLUCENT



'TI FIT VALVE WITH'



HF-VELDED SEAMS



REINFORCED WOOD  
INNER SYSTEM

## PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device.

### PRODUCT INFO

- 3 handles on each side
- Removable/changable glass fiber back side
- Non IR reflective PVC
- Inner bag with chamber system
- Carry capacity: 220 kg

Warranty period: 3 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant and Non IR reflective.
- Inner bag material: Polyester
- Granules: Polystyrene granules (flame retardant)
- Head/neck part: Birch plywood

### SPECIFICATION

- Length: 200 cm
- Width: 136 (Widest part)/74 (Head part) cm
- Height: 3,5 cm
- Weight 8,5 kg

### TEST SPECIFICATIONS

- Working temperature of the mattress is -30°C to +70°C

### ORDER INFORMATION

Art no: 271513801

RESCUE LITE FULL BODY VACUUM SPLINT MOD  
READY2GO KIT, footrest Vacuum Pump, repair kit, 3-in-1 Mattress Bag MOD

### SPARE PARTS

Art no: 12633800000

Art: Germa vacuum valve Black

Art no: 15005002817

Art: Repair Kit

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

www.germa.se  
Phone: +46 (0)44 12 30 30  
info@germa.se

## BENEFITS

## RESCUE LITE FULL BODY SPLINT MOD

The PVC is made in a Non IR reflective material, specially for the military use.

Allows preliminary X-ray & MRI diagnostics, direct in to x-ray & MRI machines with no extra handling of patient.

Glass fiber back side to improve durability and ability to slide the Full Body Splint on the ground/floor/stairs etc.

Head, neck and upper back part with wood part to secure stronger immobilisation.

Germa improved chamber technology for quick and easy preparation. Give rescuers quick handling and excellent immobilisation for the patient ensuring the granulate is evenly dispersed to give the patient stable and optimal immobilisation.

Easy to lift patient in/out from Basket Stretcher.

Light weighted and small storage, folds easily, 30% less storage space required compared to standard vacuum mattress. Easy to carry in the specially designed 3-in-1 Bag.

Robust welding, fully welded without any open seams or textiles. Making hygienic cleaning of the Full Body Splint easy.

Valve with special designed hole, ensuring almost all suction devices fit directly into the valve, no adapters are needed. Reverse air lock system, no closing of valve at stressful site.

V-shaped design to make sure upper wide part will make all sized patients the very best immobilisation of the back, neck and head.

X-designed restraint belt fixation close to the patient, support immobilisation process and comfortable for patient. Supplied with five green straps with fix lock, helping the rescuer to place the straps correctly and easy to position and adjust.

Robust and efficient material to hold vacuum for more than 24 hours, safety for patient and make confiding immobilisation.

Load capacity 220 kg, safe to handle heavy duty people.

PVC coated fabric make it easy to clean and disinfect with most common disinfectants.

Repair Kit easy to use for effective repair, quick to repair no return needed.

Deliver consistent and good immobilisation of the patient in combination with a Basket stretcher.

High quality footrest aluminium vacuum pump evacuate with more than -700 hpa, durable, quick and time saving.