



AB GERMA®

STAIR CHAIR BACK SUPPORT



PRODUCT SHEET

STAIR CHAIR BACK SUPPORT

The Stair Chair Back Support prevents the patient's head to fall through in between the handle bar when extended. The Stair Chair Back Support gives the patient a safer and more comfortable ride, making it easier for the patient to relax.

The new Stair Chair Back Support is easily attached to the chair with hooks & loops straps and it flawlessly adjusts to the handles when extended or shortened.

PRODUCT SPECIFICATIONS

PRODUCT INFO

- Mesh material.
- 6 hooks & loops straps.
- Low weight
- Easy to clean
- Give a safer transport for the patient
- The Mesh should not be folded but rolled together.
- Clean with a mild detergent and wipe off with wet cloth

Warranty period: 1 years

MATERIAL

- Mesh: PVC coated Polyester
- Straps: Polyester

SPECIFICATIONS

- Length: 57 cm
- Width: 22/13 cm
- Weight: 0,1 kg

ORDER INFORMATION

Art no: 15000000126

Stair Chair Back Support

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se