



AB GERMA®

STEP BY STEP TRANSFER  
BLACK



# PRODUCT SHEET

Evacuation tool for ergonomic and safe patient evacuation.

Intended use for patient handling on a Transfer Mattresse or Germa vacuum mattress.

When you are in a emergency situation and have to make an evacuation Step by Step straps are designed tool easing the handling of an patient. Patient handling and evacuation on flat ground, down stair cases and/or over obstacles which one must evacuate by or over.

When using Step by Step Straps for patient handling you will be evacuating in an ergonomic and safe way. A great tool that shall be used with prior training and experience by the operators.

Easy to store and to assembly on to a mattress to be used in emergency situation—slide out the patient, down stair cases (feet first). Mattress in-built friction revers brakes when sliding down in stair cases. By lifting gently at foot end you can control the speed and help to break.

## PRODUCT SPECIFICATIONS

### PRODUCT INFO

- 4 restrain strap units
- 2 steering straps with carbiner hooks and rings, one longer (head side) and one shorter (foot side)
- 1 storage bag/foot protection
- Max load capacity: 250 kg

Warranty period: 2 years

### MATERIAL

- Straps: Polyester
- Carbiner and rings: Steel

### SPECIFICATIONS

- Length: 30 cm
- Width: 20 cm
- Hight: 15 cm
- Weight: 1,6 kg

### ORDER INFORMATION

Art no: 15000000117

STEP BY STEP TRANSFER BLACKSTEP BLACK

## STEP BY STEP TRANSFER BLACK

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)  
Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)