



AB GERMA®

**STRETCHER MATTRESS
MONDIAL WITH INTEGRATED
SIDESUPPORT**



PRODUCT SHEET

The Mondial Mattress is developed in co-operation with Ferno International focusing on safety and comfort for the patient.

This mattress is designed with help of trauma teams to secure efficient handling, safety and comfort.

The mattress has built in Sensus foam for extra comfort and stability of the patient.

STRETCHER MATTRESS MONDIAL WITH INTEGRATED SIDESUPPORT



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 2 years

MATERIAL

- Outer material: Padova PVC- Polyester fibre. flame retardant
- Inner material: 70 mm V - shaped

SPECIFICATION FOR THE MATTRESS

- Length: 182 cm
- Width: 47 cm
- Height: 8 cm
- Weight 4,2 kg

ORDER INFORMATION

Art no: 271664001

Stretcher Mattress Mondial with integrated sidesupport

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

MDR 2017/745 Class 1 Medical device

Easy to assembly and adjust stretcher restraint system.

The PVC fabric is made strong and durable for long and hard use.

The Mondial mattress has been designed to fit the Ferno Mondial stretcher.

Side support via V-shaped foam in different density in the mattress to create safe feeling and comfort.

Extra pillow can be added to create added comfort for patients head and neck part.

The pillow is designed to be adjusted in position and to bridge when using head extension part of the stretcher.

Fixation of the Mondial mattress on the stretcher has two options, either with plastic hook or with Velcro straps.

For patient safety the mattress should be secured on the stretcher with the stretchers restraint system.

It is hygienic as it is fully welded and easy to clean with most common disinfectants.

STRETCHER MATTRESS MONDIAL WITH INTEGRATED SIDESUPPORT