



AB GERMA®

TRANSFER MATTRESS 12 HANDLES MOD



PRODUCT SHEET

Transfer Mattress 12 Handles MOD has been used in the ambulance service for many years with very good experience of reducing the hassle, decrease forces required in evacuation and patient handling, ergonomic and handy. The construction is very durable and strong so it can be used in ruff terrain such as a forrest and other outdoor areas and not only on smooth indoor areas.

The Transfer Mattress also has a beneficial use as a comfort/decrease pressure wounds on soldiers. It makes sure wounded soldiers who need to be transported a long distance on for example a NATO stretcher either in busses or in aeroplanes avoid pressure wounds. The mattress is thinner than a regular stretcher mattress but still gives enhanced comfort for the patient.

TRANSFER MATTRESS 12 HANDLES MOD



PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- 6 handles on each side
- Glass fiber back side

Warranty period: 2 years

MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant and Non IR reflective
- Inner material: 30 mm thick Polyurethane Foam

SPECIFICATION

- Length: 187 cm
- Width: 65 cm
- Height: 3,0 cm
- Weight 4,5 kg

ORDER INFORMATION

Art no: 271008880

TRANSFER MATTRESS 12 HANDLES MOD

SPARE PARTS

Art no: 15005002817

Art: Repair Kit

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

Made to correspondent to EN 1789

TRANSFER MATTRESS 12 HANDLES MOD

Easy to carry or slide over a body on the mattress with 12 strong handles, 6 on each long side.

Handle parts can be folded under the mattress on the stretcher.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher.

Should be used with restraints to secure the patient for safety reasons.

Comfort, ethical and ergonomic handling of the patient.

Hygienic; fully welded and easy to clean with most common disinfectants.

Fits together with Germa Step by Step evacuation Kit