




**AB GERMA®**

## TRANSFER MATTRESS 12 HANDLES



# PRODUCT SHEET

Transfer Mattress 12 Handles has been used in the ambulance service for many years with very good experience of reducing the hassle, decrease forces required in evacuation and patient handling, ergonomic and handy. The construction is very durable and strong so it can be used in ruff terrain such as a forrest and other outdoor areas and not only on smooth indoor areas.

The Transfer Mattress also has a beneficial use as a transfer mattress on a stretcher to give comfort/decrease pressure wounds. The mattress is thinner than a regular stretcher mattress but still gives enhanced comfort for the patient.

## PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- 6 handles on each side
- Glass fiber back side
- Carry capacity 250 Kg

Warranty period: 2 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- Inner material: 30 mm thick Polyurethane Foam

### SPECIFICATION

- Length: 187 cm
- Width: 65 cm
- Height: 3,0 cm
- Weight 4,5 kg

# TRANSFER MATTRESS 12 HANDLES



### ORDER INFORMATION

Art no: 271008800

TRANSFER MATTRESS 12 HANDLES Black

### SPARE PARTS

Art no: 15000000121

Art: REPAIR KIT MATTRESSES BLACK

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)  
Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)

## BENEFITS

# TRANSFER MATTRESS 12 HANDLES

Made to correspondent to EN 1789 and 1865

Easy to carry or slide over a body on the mattress with 12 strong handles, 6 on each long side.

Handles part fold under the mattress on the stretcher.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher.

Should be used with restraints to secure the body on top for safety reasons.

Comfort, ethical and ergonomic handling of the body.

Hygienic full welded easy to clean with most common disinfectants.