



AB GERMA®

TRANSFER MATTRESS  
FLEX



## PRODUCT SHEET

Transfer Mattress Flex is designed in co-operation with air ambulance professionals to secure a perfect fit for ambulance stretchers. We focused on a perfect fit in the aircraft and the ergonomics of patient handling, patient safety and comfort for patient.

## TRANSFER MATTRESS FLEX



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

Made to correspondent to EN 1789 and EN 1865  
MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- 5 handles on each side
- Glass fiber back side
- Carry capacity 250 Kg

Warranty period: 2 years

### MATERIAL

- Outer material: PVC-coated polyester fibre, flame retardant.
- PP-plate and glas fibre material to reduce friction
- Inner material: 40 mm V - shaped Polyurethane Foam + 15 mm Sensus Foam

### SPECIFICATION FOR THE TRANSFER MATTRESS FLEX

- Length: 190 cm
- Width: 65 cm
- Height: 6 cm/10 cm
- Weight 7,8 kg

### TEST SPECIFICATIONS

- Working temperature of the mattress is -30°C to +70°C

### ORDER INFORMATION

Art no: 271070000  
Transfer Mattress Flex

## BENEFITS

Made to correspondent to EN 1789  
MDR 2017/745 Class 1 Medical device

The PVC fabric is strong and durable for long and hard use.

High level of patient comfort due to the V-shaped foam plus the Sensus memory foam.

5 handles on each side, allow paramedics to easy lift, carry and slide the mattress with patient on it.

For easier patient access during transport the handle part can be folded under the mattress on the stretcher.

Sliding surface on the reverse of glas fibre material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher or from stretcher to hospital bed.

It is very easy to apply the stretcher's restraint system due to the easy access slits at the shoulder level on the Transfer Mattress Flex.

For patient safety the mattress should be secured on the stretcher with the stretcher's restraint system.

It is hygienic as it is fully welded and easy to clean with most common disinfectants.

## TRANSFER MATTRESS FLEX