



AB GERMA®

TRANSFER SHEET  
HEAVY DUTY



# PRODUCT SHEET

Heavy Duty Carry Slide sheet designed to support and help the paramedics to transfer, slide carry a patient from narrow environments and/or bed or floor on to an ambulance stretcher and in next step from the stretcher to a hospital bed. The purpose is to support handling of heavy-duty patient with safety and improved ergonomic, appreciated ergonomic tool for heavy-duty patient handling.

## TRANSFER SHEET HEAVY DUTY



FLAME RETARDANT



X-RAY MRI AND  
CT TRANSLUCENT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

MDR 2017/745 CLASS 1 MEDICAL DEVICE.

### PRODUCT INFO

- 20 handles provide optimal balance when lifting and carrying a patient onboard.

Warranty period: 2 years

### MATERIAL

- PVC: Polyester fabric material (flame retardant)
- Fiberglass reverse

### SPECIFICATIONS

- Length: 205 cm
- Width: 143 cm
- Weight 4,3 kg
- Carry capacity: 380 kg

### ORDER INFORMATION

Art no: 005000005

Transfer Sheet Heavy Duty

## BENEFITS

- 20 side handles enable safe and ergonomic lifting and carrying.
- Sliding surface on the back of the sheet, reduction of friction when sliding.
- Easy to clean fully welded product.
- Durable and strong material and design
- X-ray and MRI translucent.

### Visiting address

AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

### Post address

AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

[www.germa.se](http://www.germa.se)

Phone: +46 (0)44 12 30 30  
[info@germa.se](mailto:info@germa.se)