



**AB GERMA®**

## UNI EASY TRANSFER MATTRESS



## PRODUCT SHEET

Uni Easy Transfer mattress has been developed in co-operation with professional paramedics with focusing on hygienic, ergonomic patient handling, patient safety and comfort.

Uni Easy Transfer mattress has been tested together with medical trained trauma teams to secure efficient handling plus quick and easy transfer of patients.

The Uni Easy Transfer mattress has been designed to fit most common stretchers available on the market.

## UNI EASY TRANSFER MATTRESS



X-RAY MRI AND  
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

## PRODUCT SPECIFICATIONS

Approved according to European Standard  
EN 1789  
MDR 2017/745 Class 1 Medical device

### PRODUCT INFO

- 5 handles on each side
- Glass fiber back side
- Carry capacity 250 Kg

Warranty period: 2 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant.
- PP-plate and glass fibre material to reduce friction
- Inner material: 40 mm V - shaped Polyurethane Foam + 15 mm Sensus Foam

### SPECIFICATION FOR THE UNI EASY TRANSFER MATTRESS

- Length: 195 cm
- Width: 50 cm
- Height: 6 cm
- Weight 4,7 kg

### ORDER INFORMATION

Art no: 271030001

Uni Easy Transfer Mattress

## BENEFITS

# UNI EASY TRANSFER MATTRESS

Made to correspondent to EN 1789  
MDR 2017/745 Class 1 Medical device

The Uni Easy Transfer mattress has been designed to fit most common stretchers available on the market.

High level of patient comfort due to the V-shaped foam plus the Sensus memory foam.

To give extra comfort the PVC fabric is made from a soft and comfortable material, yet it is durable for long and hard use.

5 handles on each side, allow paramedics to easy lift, carry and slide the mattress with patient on it.

For easier patient access during transport the handle part can be folded under the mattress on the stretcher.

Sliding surface on the reverse of fiberglass material, reducing friction and allow sliding of the mattress e.g. on the floor/ground or from stretcher to stretcher or from stretcher to hospital bed.

It is very easy to apply the stretcher's restraint system due to the easy access slits at the shoulder level on the Uni Easy Transfer mattress.

For patient safety the mattress should be secured on the stretcher with the stretcher's restraint system.

It is hygienic as it is fully welded and easy to clean with most common disinfectants.