



AB GERMA®

VACUUM PILLOWS



PRODUCT SHEET

Patient positioning and fixation!

Germa Vacuum Pillows provide customized, semipermanent stabilisation, positioning and fixation of a patient in a multitude of situations, for example, dental- eye- hand and general surgery, fixation of fractured limbs, fixation at X-ray etc.

The Vacuum Pillow forms a firm support for any part of the body, the weight evenly distributed, reducing the risk of pressure points. Can be closed to accommodate structural deformities and can also be used for temporary posture correction.

PRODUCT SPECIFICATIONS

MDR 2017/745 Class 1 Medical device

PRODUCT INFO

- The Vacuum Pillows consist of an airtight shell containing granule filled inner bag.
- Available in several sizes up to mattress size
- Density and firmness can be easily adjusted
- White PVC - thinner/smooth PVC which makes it a more comfortable for the patients but still durable.
- Blue PVC - strong and durable.
- Orange PVC - same as Blue PVC

Warranty period: 2 years

MATERIAL

- Outer material: Orange: PVC-coated Polyester fibre. White and Blue PVC Foile.
flame retardant
- Inner bag: Nylone
- Granules: Polystyrene

VACUUM PILLOWS



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



HF-VELDED SEAMS

ORDER INFORMATION

White Vacum Pillows

	Art. No
Positioning Pillow 40X20cm White	22396000007
Positioning Pillow 40X30cm White	22396800001
Positioning Pillow 40X50cm White	22396600000
Positioning Pillow 55X30cm White	22393300000
Positioning Pillow 55X65cm White	22394100000
Positioning Pillow 55X85cm White	22395800000
Positioning Pillow 80X100cm White	22391100000
Positioning Pillow 80X200cm White	22392200000

Blue Vacuum Pillows

	Art. No
Positioning Pillow 40X50cm Blue	22396640000
Positioning Pillow 55X30cm Blue	22393340000
Positioning Pillow 55X65cm Blue	22394140000
Positioning Pillow 55X85cm Blue	22395840000
Positioning Pillow 80X100cm Blue	22391140000
Positioning Pillow 80X200cm Blue	22392240000

Orange Vacum Pillows

	Art. No
Positioning Pillow Rescue 55X65 Cm	22220800000
Positioning Pillow Rescue 55X85 Cm	22221600000

Vacuum Pumps

Vacuum Hand Pump	272000701
Vacuum foot pump	271000701

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se