



AB GERMA®

VACUUM SPLINT KIT AS190 MOD



PRODUCT SHEET

The Germa AS190 MOD Vacuum splints are intended to stabilise injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an optimal and comfortable immobilisation and support when evacuating the air out from the product.

All Germa MOD products are made in a Non IR reflective PVC.

VACUUM SPLINT KIT AS190 MOD



PRODUCT SPECIFICATIONS

The product is in conformity with Council MDR 2017/745 Class 1 concerning Medical Devices.

PRODUCT INFO

Wrist Splint
Arm Splint
Leg Splint
Repair kit
Hand operated Pump
Carry bag.

Warranty period: 1 years

MATERIAL

Outer material: Polyester fibre, PVC-coated, flame retardant and Non IR reflective
Inner bag: Polyester
Granules: Polystyrene granules (flame retardant)

SPECIFICATION FOR THE SPLINT SET

Length: 83 cm
Width: 32 cm
Height: 32 cm
Weight 3,85 kg

Working temperature of the mattress is -30°C to +70°C

ORDER INFORMATION

Art no: 005002810
VACUUM SPLINT KIT AS 190 MOD
– including hand operated Vacuum Pump, repair kit, bag

SPARE PARTS

Art no: 15005002817
Art: Repair Kit Olive Green

Art no: 12633800000
Art: Germa vacuum valve Black

TEST SPECIFICATIONS

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Post address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

VACUUM SPLINT KIT AS190 MOD

The hand pump is one hand operated.
Leaves you with one hand free to hold the splint and arm/leg while applying the vacuum.

High volume pump which is quick to use.
Efficient so it saves time.

Long lasting and durable material.

Multi fit valve.
Use almost any available pump.

Hooks & Loops fastener welded on the material.

The full leg splint is to ensure one size fits all patients.