



AB GERMA®

## GERMA VACUUM SPLINT SET MOD



# PRODUCT SHEET

Germa line of Vacuum Splints MOD:

The Germa Vacuum Splints MOD's are built on the original high quality, long-life models that you are used to.

Germa Vacuum Splints MOD are intended and designed to deliver premium stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

Designed and constructed to reduce needed preparation and is very reliable in its function.

All Germa MOD products are made in a Non IR reflective PVC.

## GERMA VACUUM SPLINT SET MOD



## PRODUCT SPECIFICATIONS

The product is in conformity with Council MDR 2017/745 Class 1 concerning Medical Devices.

### PRODUCT INFO

- Arm Splint
- Leg Splint
- Full Leg Splint
- Repair kit
- Hand operated Pump
- Carry bag.

Warranty period: 3 years

### MATERIAL

- Outer material: Polyester fibre, PVC-coated, flame retardant and Non IR reflective.
- Inner bag: Polyester
- Granules: Polystyrene granules (flame retardant)
- Carry Bag: Polyester

### SPECIFICATION FOR THE SPLINT SET

- Length: 83 cm
- Width: 32 cm
- Height: 32 cm
- Weight 4,75 kg

### TEST SPECIFICATIONS

- Maximum shrinkage on length: 1%
- Maximum shrinkage on width: 3%
- Working temperature of the mattress is -30°C to +70°C

### ORDER INFORMATION

Art no: 272005800

GERMA VACUUM SPLINT SET+FULL LEG VACUUM SPLINT MOD

– including hand operated Vacuum Pump, repair kit, bag

### SPARE PARTS

Art no: 15005002817

Art: Repair Kit

Art no: 12633800000

Art: Germa vacuum valve Black

Visiting address  
AB Germa  
Industrig. 54-56  
SE-291 36 Kristianstad, Sweden

Post address  
AB Germa  
Box 2127  
SE-291 02 Kristianstad, Sweden

www.germa.se  
Phone: +46 (0)44 12 30 30  
info@germa.se

## BENEFITS

## GERMA VACUUM SPLINT SET MOD

Unique inner bag with channel system, ensuring granules are evenly dispersed and quicker to use.

Robust welding, ensuring long life quality.

Leg splints with extra hooks and loops for strap allowing the splint to be used as vacuum mattress for babies and children.

Full leg splint is provided with additional hooks & loops straps to ensure one size fits all.

Hooks & loops is sewn onto a piece of PVC coated fabric to reinforce the strength of the hook & loop straps.

High volume pump which is quick to use. Efficient so it saves time.

Multi fit valve.

Use almost any available pump.