



PRODUCT SHEET

The AS Line of splints:

The AS Line is a Ferno Low-price range of splints manufactured by Germa, however, regardless of lower price they still have Germas high quality and long-life.

The AS Line splints are intended and designed to deliver stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

VACUUM SPLINTS AS LINE



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

Compliant to MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 1 years

MATERIAL

Outer material:

Polyester fibre, PVC coated (flame retardant)

Granules: Styropor granules (flame retardant)

ORDER INFORMATION:

Art: Vacuum Splint AS 140

Art no: 005002014

Immobilisation of full leg, hip or shoulder injuries to adults. Also works as a vacuum mattress for small children or babies.

Length: 108 cm

Width: 69 cm

Thickness: 4 cm

Weight: 1.32 kg



Art: Vacuum Splint AS 120

Art no: 005002013

Immobilisation of lower leg on adult and full leg, hip or shoulder on child. Also works as a vacuum mattress for babies.

Length: 74 cm

Width: 61 cm

Thickness: 4 cm

Weight: 0,90kg



Art: Vacuum Splint AS 100

Art no: 005002012

Immobilisation of arm, forearm and wrist fractures.

Length: 53 cm

Width: 35 cm

Thickness: 4 cm

Weight: 0.40 kg



Art: Vacuum Hand Pump

Art no: 272000701

One-hand-operated Vacuum Pump.

Weight: 0.6 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

MDR 2017/745 Class I Medical device

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Hooks and Loops strap allowing the splint to be used as vacuum mattress for babies and children.

Germa valve with "reverse air-lock".

X-ray compatible and translucent

VACUUM SPLINTS AS LINE

ORDER INFORMATION SETS:

Art: Vacuum Splint Kit AS190
Art no: 005002210
Short arm, long arm, long Leg,
Repair Kit, Pump and Bag
Length: 82 cm
Width: 32 cm
Thickness: 18 cm
Weight: 4.30 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se