



PRODUCT SHEET

Germa Blue Line of splints:

Blue Line is Germa's updated Low-price range of splints, however, regardless of lower price they still have Germa's high quality and long-life.

With the new Blue Line we wanted to make it as flexible as possible in respect to your needs. Therefore we have three standard sets for you to choose from to meet your needs.

The Blue Line splints are intended and designed to deliver stabilisation to injured patients during transport. They are suited for fixation of patients with hand, wrist, arm, leg, shoulder and/or hip injuries.

The splint moulds around the individual fracture giving an absolute optimal and very comfortable support to the patient.

VACUUM SPLINTS BLUE LINE



X-RAY MRI AND
CT TRANSLUCENT



FLAME RETARDANT



MULTI FIT VALVE WITH LID



HF-WELDED SEAMS

PRODUCT SPECIFICATIONS

Compliant to MDR 2017/745 Class 1 Medical device

PRODUCT INFO

Warranty period: 2 years

MATERIAL

Outer material:

Polyester fibre, PVC coated (flame retardant)

Granules: Styropor granules (flame retardant)

ORDER INFORMATION:

Art: Vacuum Splint Full Leg

Art no: 272701000

Immobilisation of full leg, hip or shoulder injuries on adults. Also works as a vacuum mattress for small children or babies.

Length: 135 cm

Width: 74 cm

Thickness: 3 cm

Weight: 1.9 kg



Art: Vacuum Splint Lower Leg

Art no: 272702000

Immobilisation of lower leg on adult and full leg, hip or shoulder on child. Also works as a vacuum mattress for babies.

Length: 96 cm

Width: 60 cm

Thickness: 3 cm

Weight: 1.2 kg



Art: Vacuum Splint Arm

Art no: 272703000

Immobilisation of arm, forearm and wrist fractures.

Length: 70 cm

Width: 40 cm

Thickness: 3 cm

Weight: 0.70 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

ORDER INFORMATION:

Art: Vacuum Splint Forearm
Art no: 272707000
Immobilisation of wrist fractures
Length: 48 cm
Width: 37 cm
Thickness: 3 cm
Weight: 0,55 kg



Art: Vacuum HandPump
Art no: 272000701
One-hand-operated
Vacuum Pump.
Weight: 0,6 kg



Art: Bag for Blue Line Splintset
Art no: 271000708
Length: 136 cm
Width: 73 cm
Thickness: 3 cm
Weight: 0,3 kg



Different sets of Splints

Art: Vacuum Splint Set Blue Line FLAF
Art no: 272708000
Full Leg, Arm, Lower Arm, Repair Kit, Pump and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 4 kg

Art: Vacuum Splint Set Blue Line LA
Art no: 272704000
Leg, Arm, Repair Kit, Pump and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 2,7 kg

Art: Vacuum Splint Set Blue Line FLA
Art no: 272705000
Full Leg, Leg, Arm, Repair Kit, Pump and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 4,6 kg

You have four different combination of kits to choose from but if a specific need please contact us.

Art: Vacuum Splint Set Blue Line LLA
Art no: 272709000
Full Leg, Leg, Arm, extension straps Repair Kit, Pump and Bag
Length: 75 cm
Width: 37 cm
Thickness: 18 cm
Weight: 4,8 kg



Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se

BENEFITS

MDR 2017/745 Class I Medical device

Robust weldings. Trustful.

PVC-coated fabric that is easy to clean and disinfect.

Full Leg splints with extra Hooks and Loops strap allowing the splint to be used as vacuum mattress for babies and children.

The Splint sets are configured after most common needs.

VACUUM SPLINTS BLUE LINE

Visiting address
AB Germa
Industrig. 54-56
SE-291 36 Kristianstad, Sweden

Postal address
AB Germa
Box 2127
SE-291 02 Kristianstad, Sweden

www.germa.se
Phone: +46 (0)44 12 30 30
info@germa.se